



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Speculum,vaginal,Graves,95x35mm

**Product Code :**

LBNY-0019-220051



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Speculum, vaginal, self-retaining.

Double blade with a screw mechanism to keep the blades open.

To dilate the vagina to examine the vaginal canal and cervix.

Material: Austenitic steel (non-quenched, non-magnetic steel).

Austenitic steel composition: 16 - 18% chromium; 10 - 14% nickel; 2 - 3% molybdenum.

Blade length: approx. 95 mm.

Blade width: approx. 35 mm.

Reusable.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

