



Email : sales@laboratoryinstrumentindia.com

Product Name :
Speculum,vaginal,Graves,95x35mm

Product Code :
LBNY-0019-220051



www.laboratoryinstrumentindia.com

Description :

Speculum, vaginal, self-retaining.
Double blade with a screw mechanism to keep the blades open.
To dilate the vagina to examine the vaginal canal and cervix.
Material: Austenitic steel (non-quenched, non-magnetic steel).
Austenitic steel composition: 16 - 18% chromium; 10 - 14% nickel; 2 - 3% molybdenum.
Blade length: approx. 95 mm.
Blade width: approx. 35 mm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

