



Email : sales@laboratoryinstrumentindia.com

Product Name :
Retractor,vaginal,38x80mm

Product Code :
LBNY-0019-2200053



www.laboratoryinstrumentindia.com

Description :

Material: Austenitic steel (non-quenched, non-magnetic steel).
To spread and open the posterior wall of the vagina to visualize the uterine cervix.
Static function instrument: vaginal speculum.
With weight which fits on the shank of the blade.
Austenitic steel composition: 18 to 20% chromium, 8 to 10% nickel.
Blade width: approx. 38 mm.
Blade length: approx. 80 mm blade.
Total length: approx. 24 cm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 40% RH to 70% RH.
Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

