



Email : sales@laboratoryinstrumentindia.com

Product Name :

Spatula,abdo,malleable,250mm

Product Code :

LBNY-0019-220054



www.laboratoryinstrumentindia.com

Description :

To hold open fatty tissue, muscle or viscera after an abdominal incision.

Static function instrument: depressor spatula.

Material: Austenitic steel (non-quenched, non-magnetic steel).

Austenitic steel composition: 16-18% chromium, 8-14% nickel, 2-3% molybdenum.

Rectangular with rounded ends.

Malleable: can be bent according to the operative field to be kept open.

Size approx. 25 x 250 mm.

Reusable.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

