



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Spatula,abdo,malleable,250mm

**Product Code :**  
LBNY-0019-220054



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

To hold open fatty tissue, muscle or viscera after an abdominal incision.  
Static function instrument: depressor spatula.  
Material: Austenitic steel (non-quenched, non-magnetic steel).  
Austenitic steel composition: 16-18% chromium, 8-14% nickel, 2-3% molybdenum.  
Rectangular with rounded ends.  
Malleable: can be bent according to the operative field to be kept open.  
Size approx. 25 x 250 mm.  
Reusable.

**Technical Specification :**

Environmental conditions:  
Humidity: 30% RH to 70% RH.  
Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

