



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Measure Tape

**Product Code :**  
LBNY-0004-820033



**Description :**

Measure Tape

**Technical Specification :**

Measure Tape

This measuring tape is ideal for DIY, domestic and industrial use. It is perfect for measuring short distances & can be operated with one hand.

It is designed for storing and carrying a variety of hand tools and power tools. Its heavy-duty construction keeps the tools secured and protected while carrying them. The box has a transparent lid on the top for easy view and access.

The multi-purpose tool box features an internal tray with handle, which enables the user to place few tools that are often needed. This easy to use tray helps in keeping the tools accessible at hand.

Material: Plastic and stainless steel

Scale: Metric and Inch 10ft (3m)

Tape thickness: 0.5mm



**Laboratory Instrument India**

