



Email : sales@laboratoryinstrumentindia.com

Product Name :
Catheter, Foley, CH10, ster, disp

Product Code :
LBNY-0021-220001



www.laboratoryinstrumentindia.com

Description :

Urethral catheter with balloon: Foley catheter 2-way with nelaton tip.
Central channel for urinary drainage, straight rounded distal end with lateral eyes, proximal end with cup connector allowing catheter to be connected to urine bag.
Diameter expressed in Charriere, French gauge.
Length expressed in cm.
Balloon capacity expressed in ml.
Side channel for inflating balloon, ending by a non-return valve with Luer tip connector.
Colour code / external diameter visible on connector.

Technical Specification :

Material: Natural latex, silicone coated.
Initial sterilisation method: ethylene oxide gas.
Balloon capacity: 3 – 5ml
Single-use.
Size selected:
Length: 20 - 30cm
Diameter: CH10
Sterile.

Weight/Volume/Dimensions
Estimated volume: 0.280cdm
Estimated weight: 0.020kg



Laboratory Instrument India