



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Medical tape,2.5cmx5m,roll

**Product Code :**  
LBNY-0021-220003



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Textile strip: woven acetate taffeta.  
Adhesive: Mixture of rubber, resins and lanolin. Traditionally incorporates zinc oxide.  
Size selected: Medical adhesive tape, roll.  
Medical adhesive tape, roll.  
High cutaneous tolerance.  
Must adhere strongly when applied to the skin, but can be removed without causing significant lesions.  
Colour: White or flesh coloured.  
Components: Aerated and / or perforated textile strip impregnated with adhesive.  
Can be torn by hand.  
Impermeable to water.  
Non-stretch.  
Length: 5m.  
Width: 2 - 2.5cm.  
Non-sterile.  
Single use.

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 0.144cdm  
Estimated weight: 0.038kg



**Laboratory Instrument India**