



Email : sales@laboratoryinstrumentindia.com

Product Name :

Medical tape,perforated,10cmx5m,roll

Product Code :

LBNY-0021-220004



www.laboratoryinstrumentindia.com

Description :

High cutaneous tolerance.

Non-stretchable.

Able to be torn by hand.

Medical adhesive tape, perforated, roll.

Impermeable to water.

Material: Textile strip: woven acetate taffeta.

Adhesive: Mixture of rubber, resins and lanolin. Traditionally incorporates zinc oxide.

Must adhere strongly when applied to the skin, but can be removed without causing any significant lesion.

Colour: White or flesh coloured.

Components: Aerated and / or perforated textile strip impregnated with adhesive, with easily detachable protective polyethylene film.

Size selected: Medical adhesive tape, roll.

Length: 5m.

Width: 10cm.

Non-sterile

Single use.

Technical Specification :

Weight/Volume

Estimated volume: 0.801cdm

Estimated weight: 0.169kg



Laboratory Instrument India