



Email : sales@laboratoryinstrumentindia.com

Product Name :
Syringe,disp,5ml,ster/BOX-100

Product Code :
LBNY-0021-220006



www.laboratoryinstrumentindia.com

Description :

Syringe, two pieces: barrel with Luer nozzle and piston/plunger or three pieces: barrel with luer nozzle, stopper and piston/plunger.
Intended use: inject fluids into, or withdraw fluids from, the body.
Blister packed.

Technical Specification :

Syringe capacity: 5ml
Syringe Material: medical grade plastic (Polyethylene (PE)/ polypropylene (PP)/ polystyrene (PS).
Syringe graduation: Graduated scale on the barrel, easy to read. Scale intervals of 0.2 or 0.5ml. Graduation lines of 1.0ml volume increment.
Barrel Nozzle: Luer lock or luer slip. The Luer nozzle shall be situated centrally.
Syringe Barrel: sufficiently transparent to allow easy measurement of the volume contained in the syringe and detection of air bubbles. Length with a maximum usable capacity of at least 10% more than the nominal capacity.
Syringe: Disposable.



Laboratory Instrument India

