



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Gown, surgical, woven

**Product Code :**  
LBNY-0021-220008



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: polyester /cotton fabric (50/50% or 65/35%) or 100% cotton fabric  
Raglan long sleeves, non-deforming cuffs in jersey: 12cm (+/- 10%)  
Washing: normal; withstands boiling and autoclaving; resists to chlorine 0.1%  
Multiple use  
Finished length of the gown: approximately: 130cm (mid-calf) (+/- 10%)  
Tie back closure and covering back panel lateral tie fastening  
Colours: blue or green  
Weight per m2: 135g-190g.

**Technical Specification :**

Gown, surgical, woven



**Laboratory Instrument India**