



Email : sales@laboratoryinstrumentindia.com

Product Name :
Gown, surgical, woven

Product Code :
LBNY-0021-220008



www.laboratoryinstrumentindia.com

Description :

Material: polyester /cotton fabric (50/50% or 65/35%) or 100% cotton fabric
Raglan long sleeves, non-deforming cuffs in jersey: 12cm (+/- 10%)
Washing: normal; withstands boiling and autoclaving; resists to chlorine 0.1%
Multiple use
Finished length of the gown: approximately: 130cm (mid-calf) (+/- 10%)
Tie back closure and covering back panel lateral tie fastening
Colours: blue or green
Weight per m2: 135g-190g.

Technical Specification :

Gown, surgical, woven



Laboratory Instrument India