



Email : sales@laboratoryinstrumentindia.com

Product Name :
Trousers,surgical,woven,size M

Product Code :
LBNY-0021-220009



www.laboratoryinstrumentindia.com

Description :

Workwear which is worn by health care worker during activities requiring a high degree of cleanliness: maternity unit, surgical unit, etc.

Technical Specification :

Material: polyester /cotton fabric (50/50% or 65/35%) or 100% cotton fabric
Trousers with drawstring waist closure
Easy fastening
Washing: normal; withstands boiling and autoclaving; resists to chlorine 0.1%
Colours: blue or green
Weight per m²: 115g-195g
Size Large (M)



Laboratory Instrument India