



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Caliber

**Product Code :**  
LBNY-0004-820034



**Description :**

Caliber

**Technical Specification :**

Caliber

Material: IP54 Waterproof Stainless Steel

With a large high-contrast LCD display, switching from mm to inch possible.

Zero switch for resetting the measured value for partial sections, very accurate and easy-to-use digital caliper.

Size: 9.5 x 3.5 x 0.6 in

Measuring Range: 0-6 inches/ 0-150mm

Accuracy:  $\pm\pm 0.03\text{mm}/0.001$

Resolution: 0.0005"/0.01mm

Battery: 3V button cell

Package Include

1\*Digital Caliper

1\* Protective Case

3\*Button Cell

1\*Screwdriver

---

# Laboratory Instrument India

