



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tunic,surgical,woven,size M

Product Code :

LBNY-0021-2200012



www.laboratoryinstrumentindia.com

Description :

Workwear which is worn by health care worker during activities requiring a high degree of cleanliness: maternity unit, surgical unit, etc.

Technical Specification :

Material: polyester /cotton fabric (50/50% or 65/35%) or 100% cotton fabric
Washing: normal; withstands boiling and autoclaving; resists to chlorine 0.1%
Multiple use
Unisex scrub top
"V-shaped" neck tunic
Easy to slip
Short sleeves
Size Medium (M)
Weight per m2: 115g-175g
Two lower patch pockets
Colours: blue or green.



Laboratory Instrument India

