



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tunic,surgical,woven, size L

Product Code :
LBNY-0021-2200013



www.laboratoryinstrumentindia.com

Description :

Workwear which is worn by health care worker during activities requiring a high degree of cleanliness: maternity unit, surgical unit, etc.

Technical Specification :

Unisex scrub top
Material: polyester /cotton fabric (50/50% or 65/35%) or 100% cotton fabric
Washing: normal; withstands boiling and autoclaving; resists to chlorine 0.1%
"V-shaped" neck tunic
Easy to slip
Two lower patch pockets
Short sleeves
Colours: blue or green
Multiple use
Size Large (L)
Weight per m2: 115g-175g



Laboratory Instrument India

