



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tunic,surgical,woven,size XL

**Product Code :**

LBNY-0021-2200014



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Workwear which is worn by health care worker during activities requiring a high degree of cleanliness: maternity unit, surgical unit, etc.

**Technical Specification :**

Material: polyester /cotton fabric (50/50% or 65/35%) or 100% cotton fabric  
Washing: normal; withstands boiling and autoclaving; resists to chlorine 0.1%  
Multiple use  
Unisex scrub top  
"V-shaped" neck tunic  
Two lower patch pockets  
Colours: blue or green  
Easy to slip  
Short sleeves  
Weight per m2: 115g-175g  
Size Extra Large (XL)



**Laboratory Instrument India**

