



Email : sales@laboratoryinstrumentindia.com

Product Name :

Coat, medical, woven, white, medium size

Product Code :

LBNY-0021-2200016



www.laboratoryinstrumentindia.com

Description :

Material: polyester /cotton fabric (50/50%, 65/35%, 67%/33%) or 100% cotton fabric
Snap closure or button front closure
Unisex medical coat with long sleeves
Two patch pockets and a chest pocket
Notched lapel collar
Washing: normal; resists to chlorine 0.1%
Mid-thigh length
Colour: white
Multiple use
Weight per m2: 175g-240g

Technical Specification :

Size: Medium (M)



Laboratory Instrument India

