



Email : sales@laboratoryinstrumentindia.com

Product Name :

Coat, medical, woven, white, large size

Product Code :

LBNY-0021-2200017



www.laboratoryinstrumentindia.com

Description :

Material: polyester /cotton fabric (50/50%, 65/35%, 67%/33%) or 100%cotton fabric

Unisex medical coat with long sleeves

Washing: normal; resists to chlorine 0.1%

Snap closure or button front closure

Two patch pockets and a chest pocket

Notched lapel collar

Multiple use

Mid-thigh length

Color: white

Weight per m2: 175g-240g

Technical Specification :

Size: Large (L)



Laboratory Instrument India

