



Email : sales@laboratoryinstrumentindia.com

Product Name :
Catheter,urethral,CH12,ster,disp

Product Code :
LBNY-0021-2200019



www.laboratoryinstrumentindia.com

Description :

Atraumatic, soft, rounded closed distal end with 2 lateral eyes, proximal end with cup connector allowing catheter to be connected with urine bag.
Urethral catheter without balloon, nelaton tip.
Length expressed in cm.
Colour code cup connector.

Technical Specification :

Material: Polyvinyl chloride (PVC).
Initial sterilization method: Ethylene oxide gas.
Single-use
Sterile.
Size selected:
Diameter CH12,
Length: 30 - 40cm.



Laboratory Instrument India

