



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vernier Calipers

Product Code :
LBNY-0004-820035



Description :

Vernier Calipers

Technical Specification :

Vernier Calipers

Material production: Made of real stainless steel hardened Chromeplated, fine ground and polished.

Accuracy? Professional precision measuring tool for industrial use (particularly on metal workpieces), totally different from DIY hand measuring one (like any carbon steel caliper, which is not good to measure metal workpieces because of the low hardness)

4 ways to use? inside/outside/step/depth measurement, which almost meets all basic precision measuring needs

Multi Function: The metric caliper is 0-6 inch / 150mm, it can measure outside diameter, inside diameter, depth and step value of an exquisite item, a multifunctional tool for you. Measuring range: 0 - 150mm (0 - 6 inch)

Slide is smooth, no grittiness. Integrated depth rod. Highly polished surface, Heavyweight, Rugged Design with Hardened Stainless Measuring Surfaces

Measuring Range: 0-6 mm

Material : Stainless steel

Accuracy: 1mm

Laboratory Instrument India

