



Email : sales@laboratoryinstrumentindia.com

Product Name :

Catheter,urethral,CH14,ster,disp

Product Code :

LBNY-0021-2200020



www.laboratoryinstrumentindia.com

Description :

Atraumatic, soft, rounded closed distal end with 2 lateral eyes, proximal end with cup connector allowing catheter to be connected with urine bag.

Urethral catheter without balloon, nelaton tip.

Colour code cup connector.

Length expressed in cm.

Technical Specification :

Initial sterilization method: Ethylene oxide gas

Material: Polyvinyl chloride (PVC).

Single-use

Sterile.

Size selected:

Length: 30 - 40cm,

Diameter CH14.



Laboratory Instrument India

