



Email : sales@laboratoryinstrumentindia.com

Product Name :
Catheter, Foley, CH12, ster, disp

Product Code :
LBNY-0021-2200021



www.laboratoryinstrumentindia.com

Description :

Urethral catheter with balloon: Foley catheter 2-way with nelaton tip.
Side channel for inflating balloon, ending by a non-return valve with Luer tip connector.
Central channel for urinary drainage, straight rounded distal end with lateral eyes, proximal end with cup connector allowing catheter to be connected to urine bag.
Balloon capacity expressed in ml.
Length expressed in cm.
Colour code / external diameter visible on connector.

Technical Specification :

Initial sterilisation method: ethylene oxide gas.
Material: Natural latex, silicone coated.
Sterile.
Single-use.
Size selected:
Length: 30 - 40cm,
Diameter CH12,
Balloon capacity: 5 – 15ml.

Laboratory Instrument India

