



Email : sales@laboratoryinstrumentindia.com

Product Name :
Catheter, Foley, CH14, ster, disp

Product Code :
LBNY-0021-2200022



www.laboratoryinstrumentindia.com

Description :

Central channel for urinary drainage, straight rounded distal end with lateral eyes, proximal end with cup connector allowing catheter to be connected to urine bag.

Urethral catheter with balloon: Foley catheter 2-way with nelaton tip.

Side channel for inflating balloon, ending by a non-return valve with Luer tip connector.

Colour code / external diameter visible on connector.

Length expressed in cm.

Balloon capacity expressed in ml.

Technical Specification :

Material: Natural latex, silicone coated.

Sterile.

Single-use.

Initial sterilisation method: ethylene oxide gas.

Size selected:

Diameter CH14,

Length: 30 - 40cm,

Balloon capacity: 5 – 15ml.

Laboratory Instrument India

