



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Catheter, Foley, CH14, ster, disp

**Product Code :**  
LBNY-0021-2200022



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Central channel for urinary drainage, straight rounded distal end with lateral eyes, proximal end with cup connector allowing catheter to be connected to urine bag.  
Urethral catheter with balloon: Foley catheter 2-way with nelaton tip.  
Side channel for inflating balloon, ending by a non-return valve with Luer tip connector.  
Colour code / external diameter visible on connector.  
Length expressed in cm.  
Balloon capacity expressed in ml.

**Technical Specification :**

Material: Natural latex, silicone coated.  
Sterile.  
Single-use.  
Initial sterilisation method: ethylene oxide gas.  
Size selected:  
Diameter CH14,  
Length: 30 - 40cm,  
Balloon capacity: 5 – 15ml.

---

# Laboratory Instrument India

