



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Simple Protractor

**Product Code :**  
LBNY-0004-820036



**Description :**

Simple Protractor

**Technical Specification :**

Simple Protractor  
Used to draw radial lines, setting bevels and transferring angles  
Made of polished stainless steel  
Double graduation from 0-180 degree  
Brass locking nut  
Scale range;0-10cm  
Angle : 0o- 180o  
Material: stainless steel



**Laboratory Instrument India**