



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bag,urine,collecting,2000ml

Product Code :
LBNY-0021-2200035



www.laboratoryinstrumentindia.com

Description :

Urine collection bag.
With reinforced eyelets for hanging;
Kink resistant and transparent plastic inlet tube, length 85-95 cm, with universal connector and protective cap.
Medical grade plastic bag;
Capacity 2000 ml;
With graduations, every 100 ml to allow proper reading of the liquid contained in the bag.

Technical Specification :

Material:
Tube & connector/protective cap: polyvinyl chloride (PVC)& Single-use
Bag: polyvinyl chloride (PVC) or ethylene vinyl acetate (EVA).
Non sterile



Laboratory Instrument India