



Email : sales@laboratoryinstrumentindia.com

Product Name :

Syringe,feeding,50ml,catheter tip,ster

Product Code :

LBNY-0021-2200037



www.laboratoryinstrumentindia.com

Description :

Clear syringe barrel, precise graduation in ml for the volume dosage, permanent marking, easy readability

Syringe, 3-pieces: barrel, plunger and gasket.

Capacity: 50ml or 60ml

Disposable

Sterile

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

Centric or eccentric catheter tip

Safe plunger backstop.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.300cdm

Estimated weight: 0.040kg



Laboratory Instrument India

