



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Syringe,feeding,50ml,catheter tip,ster

**Product Code :**

LBNY-0021-2200037



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Clear syringe barrel, precise graduation in ml for the volume dosage, permanent marking, easy readability

Syringe, 3-pieces: barrel, plunger and gasket.

Capacity: 50ml or 60ml

Disposable

Sterile

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

Centric or eccentric catheter tip

Safe plunger backstop.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.300cdm

Estimated weight: 0.040kg



**Laboratory Instrument India**

