



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube,asp/feed,CH16,L125cm,ster,disp

**Product Code :**

LBNY-0021-2200039



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Sterile enteral feeding/aspiration catheters for single use  
Proximal end with cup connector, conical tip, allowing the tube to be connected to devices such as feeding syringes (catheter tip) or suction pump system  
Single channel tube  
Open distal end with four side windows  
Colour code/ external diameter: Visible on cup connector  
Material: Polyvinyl chloride (PVC)  
Tube with markings at 40, 50, 60 and 70cm from the distal end  
Initial sterilisation method: Ethylene oxide gas  
Disposable  
Sterile  
Length expressed in cm  
Size selected: Diameter CH16, length approximately 125cm.

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 0.300cdm  
Estimated weight: 0.031kg

---

# Laboratory Instrument India

