



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,asp/feed,CH08,L125cm,ster,disp

**Product Code :**

LBNY-0021-2200041



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Sterile enteral feeding/aspirating catheter for single use.

Single channel tube.

Open distal end with four side windows.

Tube with markings at 40, 50, 60 and 70cm from the distal end.

Proximal end with cup connector, conical tip, allowing the tube to be connected to devices such as feeding syringes (catheter tip) or suction pump system.

Material: Polyvinyl chloride (PVC)

Size selected: Diameter CH08, length approximately 125cm

Disposable

Length expressed in cm.

Colour code/ external diameter: Visible on cup connector.

Sterile

Initial sterilisation method: Ethylene oxide gas.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.232cdm

Estimated weight: 0.017kg

---

# Laboratory Instrument India

