



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,asp/feed,CH12,L125cm,ster,disp

**Product Code :**

LBNY-0021-2200042



www.laboratoryinstrumentindia.com

**Description :**

Sterile enteral feeding/aspirating catheter for single use.

Proximal end with cup connector, conical tip, allowing the tube to be connected to devices such as feeding syringes (catheter tip) or suction pump system.

Single channel tube.

Colour code/ external diameter: Visible on cup connector.

Open distal end with four side windows.

Tube with markings at 40, 50, 60 and 70cm from the distal end.

Length expressed in cm.

Material: Polyvinyl chloride (PVC).

Sterile

Disposable

Initial sterilisation method: Ethylene oxide gas

Size selected: Diameter CH12, length approximately 125cm.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.232cdm

Estimated weight: 0.022kg.

---

# Laboratory Instrument India

