



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube,feeding,CH08,L40cm,ster,disp

**Product Code :**

LBNY-0021-2200043



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Sterile enteral feeding catheter for single use.

Proximal end with connector, Luer tip + stopper, allowing the tube to be connected to devices such as feeding syringes (Luer tip).

Single channel tube.

Rounded, smooth distal end with two side windows.

Material: Polyvinyl chloride (PVC).

Disposable

Sterile

Size selected: Diameter CH08, length approximately: 40 - 50cm

Tube with marking at 20cm from the distal end.

Length expressed in cm.

Colour code/ external diameter: Visible on connector.

Initial sterilisation method: Ethylene oxide gas.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.101cdm

Estimated weight: 0.005kg

---

# Laboratory Instrument India

