



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,feeding,CH05,L40cm,ster,disp

**Product Code :**

LBNY-0021-2200044



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Sterile enteral feeding catheter for single use.

Single channel tube.

Rounded, smooth distal end with two side windows.

Tube with marking at 20cm from the distal end.

Proximal end with connector, Luer tip + stopper, allowing the tube to be connected to devices such as feeding syringes (Luer tip).

Material: Polyvinyl chloride (PVC)

Size selected: Diameter CH05, length approximately 40 - 50cm

Sterile

Disposable

Length expressed in cm.

Colour code/ external diameter: Visible on connector.

Initial sterilisation method: Ethylene oxide gas.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.099cdm

Estimated weight: 0.005kg

---

# Laboratory Instrument India

