



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube,suction,CH08,L50cm,ster,disp

Product Code :

LBNY-0021-220045



www.laboratoryinstrumentindia.com

Description :

Sterile suction catheter for single use

Proximal end with cup connector, conical tip, allowing the tube to be connected to devices such as suction pump system

Open distal end, straight, with two side windows

Single channel tube

Material: Polyvinyl chloride (PVC)

Length expressed in cm

Sterile

Disposable

Size selected: Diameter CH08, length approximately 50cm

Initial sterilisation method: Ethylene oxide gas.

Colour code/ external diameter: Visible on cup connector.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.090cdm

Estimated weight: 0.040kg

Laboratory Instrument India

