



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,suction,CH10,L50cm,ster,disp

**Product Code :**

LBNY-0021-2200046



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Sterile suction catheter for single use.

Proximal end with cup connector, conical tip, allowing the tube to be connected to devices such as suction pump system.

Single channel tube.

Length expressed in cm.

Colour code/ external diameter: Visible on cup connector.

Open distal end, straight, with two side windows.

**Technical Specification :**

Material: Polyvinyl chloride (PVC).

Sterile

Disposable

Initial sterilisation method: Ethylene oxide gas

Size selected: Diameter CH10, length approximately 50cm.



**Laboratory Instrument India**

