



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,suction,CH14,L50cm,ster,disp

**Product Code :**

LBNY-0021-2200047



www.laboratoryinstrumentindia.com

**Description :**

Sterile suction catheter for single use.

Single channel tube.

Open distal end, straight, with two side windows.

Proximal end with cup connector, conical tip, allowing the tube to be connected to devices such as suction pump system.

Material: Polyvinyl chloride (PVC).

Size selected: Diameter CH14, length approximately 50cm

Sterile

Disposable

Initial sterilisation method: Ethylene oxide gas

Length expressed in cm.

Colour code/ external diameter: Visible on cup connector.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.090cdm

Estimated weight: 0.040kg

---

# Laboratory Instrument India

