



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube,suction,CH16,L50cm,ster,disp

Product Code :

LBNY-0021-220048



www.laboratoryinstrumentindia.com

Description :

Sterile suction catheter for single use.

Proximal end with cup connector, conical tip, allowing the tube to be connected to devices such as suction pump system.

Single channel tube.

Material: Polyvinyl chloride (PVC).

Size selected: Diameter CH16, length approximately 50cm

Open distal end, straight, with two side windows.

Sterile

Disposable

Initial sterilisation method: Ethylene oxide gas

Length expressed in cm.

Colour code/ external diameter: Visible on cup connector.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.090cdm

Estimated weight: 0.040kg

Laboratory Instrument India

