



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Wrench

Product Code :
LBNY-0004-820038



Description :

Pipe Wrench

Technical Specification :

Pipe Wrench

Sturdy ductile-iron housing and I-beam handle for ultimate strength and durability.

The full floating forged hook jaw allows for superior gripping and quick adjustments.

Plumbing wrench features self-cleaning threads with easily replaceable hook jaw, heel jaw and spring assembly.

Length : 10 inch

Jaw Type : Full-floating forged hook jaw with self-cleaning threads



Laboratory Instrument India