



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,endotrach,3.5,w/o cuff,ster,disp

**Product Code :**

LBNY-0021-2200050



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Endotracheal tube shall meet the required standards describing the dimensions, the tube characteristics, tube point and the tube markings.

Endotracheal tube with:

Radio opaque mark.

With Murphy's eye.

Graduations.

Open distal end and Magill-type point with oral angle of 37.5°.

Standard connector (ext. Ø 15mm) at the proximal end enabling the tube to be connected to the ventilation system.

Endotracheal tube without cuff.

Material: Polyvinyl chloride (PVC).

Disposable.

Sterile.

Size: Ø internal 3.5mm.

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

**Technical Specification :**

Tube, endotracheal, 3.5, without cuff, sterile, disposable

---

# Laboratory Instrument India

