



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,endotrach,6.5,w/cuff,ster,disp

**Product Code :**

LBNY-0021-220051



www.laboratoryinstrumentindia.com

**Description :**

Endotracheal tube shall meet the standards concerning the dimension, the tube, the point, the cuff and the marking.

Endotracheal tube with:

Cuff located near the distal end: Low pressure cuff, inflating system with valve and Luer tip connector.

Open distal end and Magill-type point with oral angle of 37.5°.

Radio opaque mark.

With Murphy's Eye.

Graduations.

Standard connector (ext. Ø 15mm) at the proximal end enabling the tube to be connected to the ventilation system.

Material: Polyvinyl chloride (PVC).

Disposable.

Sterile.

Size: Ø internal 6.5mm.

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

**Technical Specification :**

Tube, endotracheal, 6.5, with cuff, sterile, disposable

---

# Laboratory Instrument India

