



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,endotrach,7,w/cuff,ster,disp

**Product Code :**

LBNY-0021-220052



www.laboratoryinstrumentindia.com

**Description :**

Endotracheal tube shall meet the standards concerning the dimension, the tube, the point, the cuff and the marking.

Endotracheal tube with:

Standard connector (ext. Ø 15mm) at the proximal end enabling the tube to be connected to the ventilation system.

Open distal end and Magill-type point with oral angle of 37.5°.

Cuff located near the distal end: Low pressure cuff, inflating system with valve and Luer tip connector.

Radio opaque mark.

With Murphy's eye.

Material: Polyvinyl chloride (PVC).

Sterile.

Disposable.

Graduations.

Size: Ø internal 7mm.

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

**Technical Specification :**

Tube, endotracheal, 7, with cuff, sterile, disposable

---

# Laboratory Instrument India

