



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube, endotrach, 7.5, w/cuff, ster, disp

**Product Code :**

LBNY-0021-2200053



www.laboratoryinstrumentindia.com

**Description :**

Endotracheal tube shall meet the standards concerning the dimension, the tube, the point, the cuff and the marking.

Open distal end and Magill-type point with oral angle of 37.5°.

Cuff located near the distal end: Low pressure cuff, inflating system with valve and Luer tip connector.

Endotracheal tube with:

Radio opaque mark.

With Murhpy's eye.

Graduations.

Standard connector (ext. Ø 15mm) at the proximal end enabling the tube to be connected to the ventilation system.

Material: Polyvinyl chloride (PVC).

Sterile.

Disposable.

Size: Ø internal 7.5mm.

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

**Technical Specification :**

Tube, endotracheal, 7.5, with cuff, sterile, disposable

---

# Laboratory Instrument India

