



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Tourniquet, latex-free, 50cm

**Product Code :**  
LBNY-0021-2200055



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It is typically used by emergency medical services (EMS) for emergency blood loss prevention or for routine blood sampling. A strap or tubing intended to be wrapped around a patient's limb (arm or leg) and manually tightened to reduce blood circulation to or from the portion of the limb distal to where the device is applied.

**Technical Specification :**

Tourniquet used to compress a limb to stop blood flow.

Material: Latex free.

Solid elastic strip.

Reusable, non-sterile

Autoclavable at 121°C.

Length: 50-75 cm

Width: for solid design: 1.8-2.3 cm.



**Laboratory Instrument India**