



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tourniquet, latex-free, 50cm

Product Code :
LBNY-0021-2200055



www.laboratoryinstrumentindia.com

Description :

It is typically used by emergency medical services (EMS) for emergency blood loss prevention or for routine blood sampling. A strap or tubing intended to be wrapped around a patient's limb (arm or leg) and manually tightened to reduce blood circulation to or from the portion of the limb distal to where the device is applied.

Technical Specification :

Tourniquet used to compress a limb to stop blood flow.

Material: Latex free.

Solid elastic strip.

Reusable, non-sterile

Autoclavable at 121°C.

Length: 50-75 cm

Width: for solid design: 1.8-2.3 cm.



Laboratory Instrument India