



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Threading Machine

Product Code :
LBNY-0004-820039



Description :

Pipe Threading Machine

Technical Specification :

Pipe Threading Machine

Ideal for fast, clean pipe cutting by hand or power

Features an extra-long shank that protects adjustment threads

Extra-large handle allows for quick and easy adjustments

Can be converted to 3-wheel cutter by replacing rollers with cutter wheel

Pipe cutter has a 1/8-inch (3 mm) to 2-inch (50mm) pipe capacity

Converted to 3-wheel cutter by replaciug rollers with cutter wheel



Laboratory Instrument India