



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bandage,elastic,7.5cmx5m,roll

Product Code :
LBNY-0021-2200061



www.laboratoryinstrumentindia.com

Description :

Stretch bandage roll with clips.
Non-adhesive, unbleached colour and non detectable by X-ray.
Thread count:
Material: Crepe bandage produced by combining high twist cotton threads with normal twist cotton threads in warp. 100% cotton.
Weft: 54 threads/dm + or - 2 threads unbleached cotton.
Warp: 120 threads/dm, unbleached cotton, high-twist.
Size selected: Bandage elastic, roll
Length: Not stretched: +/- 3m. Stretched: +/- 5m (elasticity +/- 150 %).
Reusable.
Non-sterile.
Width: 7.5 – 8 cm

Technical Specification :

Weight/Volume
Estimated volume: 0.283cdm
Estimated weight: 0.036kg.

Laboratory Instrument India

