



Email : sales@laboratoryinstrumentindia.com

Product Name :

Cotton wool,500g,roll,non-ster

Product Code :

LBNY-0021-2200064



www.laboratoryinstrumentindia.com

Description :

High absorbent surgical grade cotton wool.
Material: 100% cotton surgical grade, purified, bleached, and carded.
Not pre-cut.
Roll presentation.
Water absorbent.
Single use.
Non-sterile.
Size selected: Cotton wool, roll of 500g.
Net weight: 500g.

Technical Specification :

Weight/Volume
Estimated volume: 3.774cdm
Estimated weight: 0.599kg



Laboratory Instrument India

