



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Gauze,roll,90cmx100m,non-ster

**Product Code :**

LBNY-0021-2200065



www.laboratoryinstrumentindia.com

**Description :**

Non sterile gauze roll with selvedge

Material: Absorbent gauze, 100% cotton.

Type of gauze: 17 threads/cm<sup>2</sup> (grammage 23g /m<sup>2</sup>).

Non-detectable by X-ray

Components: Bleached, purified textile, plain weave.

Size selected: Gauze roll (gauze folded in zigzag could be an alternative).

Length:100m.

Width: 65 to 90cm.

Non-sterile.

Single use

Thread count:

Weft: 66 to 74 threads/dm.

Warp: 95 to 105 threads/dm.

Weight: +/- 23g/m<sup>2</sup>.

**Technical Specification :**

Weight/Volume

Estimated volume: 12.996cdm

Estimated weight: 2.322kg



**Laboratory Instrument India**