



Email : sales@laboratoryinstrumentindia.com

Product Name :

Compress, gauze, 10x10cm, n/ster/PAC-100

Product Code :

LBNY-0021-2200066



www.laboratoryinstrumentindia.com

Description :

Non sterile gauze compress, folded edges, pack of 100 units.

Material: Absorbent gauze, 100 % cotton.

Non-detectable by X-ray.

Components: Bleached, purified textile, plain weave.

Type of gauze: 17 threads/cm² (grammage 23g/m²).

Number of folds (thickness): 12 ply

Non-sterile.

Single use

Size selected: Compress gauze folded.

Length: 10cm.

Width: 10cm.

Thread count:

Weft: 66 to 74 threads/dm.

Warp: 95 to 105 threads/dm.

Weight: +/- 23g/m².

Technical Specification :

Weight/Volume

Estimated volume: 2.267cdm

Estimated weight: 0.311kg



Laboratory Instrument India